

Features

- Mid-wheel drive allows for superior stability and provides a tight 18" turning radius, ideal for indoor maneuverability
- Durable construction with easy-to-access components
- Removable flip-up armrests have adjustable height, width and angle
- Programmable Dynamic controller with simple, common connectors
- Adjustable length controller mount can be affixed to either the left or right armrest and can be fitted with an optional swing-away arm
- Standard seat belt and seat back storage pouch
- Flip-up footplate with adjustable height and angle
- Elevating leg rests available (sold separately)
- Comfortable swivel Captain's seat with fold-down backrest and adjustable height headrest available in 18" & 20" widths
- All seats are height adjustable and feature 25 degrees of recline
- Batteries included

Color Options

Metallic Red

Blue



Accessories

- Elevating Leg Rests
- Swingaway Controller Arm
- Crutch/Cane Holder
- Cylinder Tank Carrier
- Armrest Bag
- Backpack
- Cup Holder

Specifications

Max Weight Capacity	300 lbs.	Power	
Max Speed	4 mph	Motor	Two 24V x 320W x 4600rpm
Max Range	20 miles	Controller	Dynamic A50A
Max Grade/Incline	12°	Batteries	Two DC12V x 35AH (UI)
Turning Radius	18"	Battery Charger	5 Amp, Offboard
Ground Clearance	3"	Brakes	Electromagnetic
Dimensions		Freewheel Mode	Yes
Overall Length	43"	Tires (Flat-free)	
Overall Width	24"	Front & Rear Casters	6" x 2"
Floor to Seat Height	20"-23"	Mid Drive Wheels	10" x 3"
Seat Dimensions		Weights	
Available Sizes	18"(w) x 18"(d) 20"(w) x 18"(d)	Without Batteries	178 lbs.
Back Height	21"	With Batteries	226 lbs.
Armrest to Seat Height	7"		

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

ActiveCare
Medical

2 Harbison Way
Columbia, SC 29212
866.621.4933
activecaremed.com

Information contained herein is correct at the time of publication. ActiveCare reserves the right to alter specifications without notice.