

Features

- Full-size power wheelchair performance in a lightweight, easily transportable package
- Can be disassembled quickly simply by pulling up on a quick release lever
- Removable flip-up armrests have adjustable width and angle
- Programmable controller with simple, common connectors
- Adjustable length controller mount can be affixed to either the left or right armrest and can be fitted with an optional swing-away arm
- Standard anti-tip wheels for added safety
- Flip-up footplate with adjustable height and angle
- Comfortable, swivel Captain's seat with fold-down backrest
- Armrest reflectors for enhanced visibility
- Batteries included

Color Options

Metallic Red

Midnight Blue



Accessories

- Elevating Leg Rests
- Backpack
- Swingaway Controller Arm
- Cup Holder
- Armrest Bag



The Cobalt disassembles easily into three pieces.

Specifications

Max. Weight Capacity 250 lbs.

Max. Range 7 miles

Max. Grade/Incline 6°

Turning Radius 29"

Ground Clearance 2"

Dimensions

Overall Length 38.5"

Overall Width 24"

Floor to Seat Height 20"-23"

Seat Dimensions

Available Size 16"(w) x 15"(d)

Back Height 15"

Armrest to Seat Height 7"

Power

Motor Two 24V x 120W x 4700rpm

Controller Dynamic A50A

Batteries Two DC12V x 21AH

Battery Charger 2 Amp, Offboard

Brakes Electromagnetic

Freewheel Mode Yes

Tires (Flat-free)

Front Casters 8" x 2"

Rear Drive Wheels 9" x 3"

Weights

Without Batteries 93 lbs.

With Batteries 121 lbs.

Heaviest pc. wt. (front section) 52 lbs.

ActiveCare
Medical

2 Harbison Way
Columbia, SC 29212
866.621.4933
activecaremed.com

Information contained herein is correct at the time of publication. ActiveCare reserves the right to alter specifications without notice.