

Features

- Full-size power wheelchair performance in a lightweight, easily transportable package
- Can be disassembled quickly and easily with the simple pull of a lever
- Removable flip-up armrests have adjustable height, width and angle
- Programmable controller with simple, common connectors
- Adjustable length controller mount can be affixed to either the left or right armrest and can be fitted with an optional swing-away arm
- Standard anti-tip wheels for added safety
- Standard seat belt and seat back storage pouch
- Flip-up footplate with adjustable height and angle
- Comfortable, swivel Captain's seat with fold-down backrest
- Armrest reflectors for enhanced visibility
- Batteries included

Color Options

Metallic Red

Midnight Blue

Accessories

- Elevating Leg Rests
- Swingaway Controller Arm
- Crutch/Cane Holder
- Cylinder Tank Carrier
- Armrest Bag
- Backpack



The Cobalt quickly and easily disassembles into 3 pieces for transport

Specifications

Max. Weight Capacity	250 lbs.	Power	Motor Two 24V x 120W x 4700rpm
Max. Range	7 miles	Controller	Dynamic A50A
Max. Grade/Incline	6°	Batteries	Two DC12V x 21AH
Turning Radius	29"	Battery Charger	2 Amp, Offboard
Ground Clearance	2"	Brakes	Electromagnetic
Dimensions		Freewheel Mode	Yes
Overall Length	38.5"	Tires (Flat-free)	
Overall Width	24"	Front Casters	8" x 2"
Floor to Seat Height	20"-23"	Rear Drive Wheels	9" x 3"
Seat Dimensions		Weights	
Available Size	16"(w) x 15"(d)	Without Batteries	93 lbs.
Back Height	15"	With Batteries	121 lbs.
Armrest to Seat Height	7"	Heaviest pc. wt. (front section)	52 lbs.

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

ActiveCare
Medical

2 Harbison Way
Columbia, SC 29212
866.621.4933
activecaremed.com

Information contained herein is correct at the time of publication. ActiveCare reserves the right to alter specifications without notice.