

Osprey™ Large Mobility Scooter

Features



- **High ground clearance and pneumatic tires ideal for outdoor travel**
- Industry-leading 500 lb. weight capacity
- 38 mile maximum range
- Adjustable speed up to a maximum of 9 mph
- Adjustable angle, delta tiller for comfortable, ergonomic control
- Standard safety package includes headlight, horn, turn signals and dual rearview mirrors
- Front and rear suspension for a more comfortable ride
- Anti-tip wheels for greater safety
- Comfortable swivel Captain's seat with fold-down backrest and adjustable height headrest available in 20" and 22" widths
- Removable flip-up armrests have adjustable height, width and angle
- All seats are height adjustable and feature 25 degrees of recline
- Standard metal wire basket
- Batteries included



Color

Metallic Gray

Accessories

- Scooter Trailer
- Crutch/Cane Holder
- Cylinder Tank Carrier
- Armrest Bag
- Backpack
- Rear Basket
- Dust Cover
- Cup Holder

Specifications

Max Weight Capacity	500 lbs.	Power	
Max Speed	9 mph	Motor	24V x 920W x 3700rpm
Max Range	38 miles	Controller	Dynamic 160A
Max Grade/Incline	12°	Batteries	Two DC12V x 75AH
Turning Radius	73"	Battery Charger	5 Amp, Offboard
Ground Clearance	7"	Brakes	Electromagnetic/Drum
Dimensions		Freewheel Mode	Yes
Overall Length	54"	Tires (Pneumatic)	
Overall Width	29"	Front & Rear	13" x 5"
Floor to Seat Height	25"-28"	Weights	
Seat Dimensions		Without Batteries	235 lbs.
Available Sizes	20"(w) x 19"(d) 22"(w) x 20"(d)	With Batteries	300 lbs.
Back Height	21"	<i>Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.</i>	
Armrest to Seat Height	7"		

ActiveCare
Medical

2 Harbison Way
Columbia, SC 29212
866.621.4933
activecaremed.com

Information contained herein is correct at the time of publication. ActiveCare reserves the right to alter specifications without notice.