

## Features

- Durable construction with easy-to-access components
- Removable flip-up armrests have adjustable height, width and angle
- Programmable Dynamic controller with simple, common connectors
- Adjustable length controller mount can be affixed to either the left or right armrest and can be fitted with an optional swing-away arm
- Standard anti-tip wheels for added safety
- Standard seat belt and seat back storage pouch
- Flip-up footplate with adjustable height and angle
- Elevating leg rests available (sold separately)
- Comfortable swivel Captain's seat with fold-down backrest and adjustable height headrest available in 18", 20" and 22" widths
- All seats are height adjustable and feature 25 degrees of recline
- Includes handy rear trunk and rear safety reflectors
- Batteries included

## Color Options

Metallic Red




Midnight Blue

## Accessories

- Elevating Leg Rests
- Swingaway Controller Arm
- Crutch/Cane Holder
- Cylinder Tank Carrier
- Armrest Bag
- Backpack
- Cup Holder

## Specifications

<b>Max Weight Capacity</b>	<b>350 lbs.</b>	<b>Power</b>	
<b>Max Speed</b>	4 mph	Motor	Two 24V x 320W x 4600rpm
<b>Max Range</b>	20 miles	Controller	Dynamic A50A
<b>Max Grade/Incline</b>	12°	Batteries	Two DC12V x 35AH (UI)
<b>Turning Radius</b>	31"	Battery Charger	5 Amp, Offboard
<b>Ground Clearance</b>	3"	<b>Brakes</b>	Electromagnetic
<b>Dimensions</b>		<b>Freewheel Mode</b>	Yes
Overall Length	37"	<b>Tires (Flat-free)</b>	
Overall Width	24"	Front Casters	8" x 2"
Floor to Seat Height	20"-23"	Rear Drive Wheels	10" x 3"
<b>Seat Dimensions</b>		<b>Weights</b>	
Available Sizes	18"(w) x 18"(d) 20"(w) x 19"(d) 22"(w) x 20"(d)	Without Batteries	143 lbs.
Back Height	21"	With Batteries	191 lbs.
Armrest to Seat Height	7"	<i>Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.</i>	



**ActiveCare**  
Medical

2 Harbison Way  
Columbia, SC 29212  
866.621.4933  
activecaremed.com

Information contained herein is correct at the time of publication. ActiveCare reserves the right to alter specifications without notice.