

Features

- Lightweight steel frame can be folded to a compact 13" width for easy portability, transport and storage
- Adjustable height, flip-up armrests
- Programmable Dynamic controller with common connectors
- Adjustable length controller mount can be affixed to either the left or right armrest and can be fitted with an optional swing-away arm
- Heavy-duty, flat-free, non-marking tires
- Standard anti-tip wheels for added safety
- Standard seat belt and seat back storage pouch
- Swing-away footrests, aluminum footplates and calf strap
- Elevating leg rests available (sold separately)
- Padded, fold-down backrest and cushioned seat available in 20" and 22" widths
- Manual wheel locks for additional safety

Specifications

Maximum Weight Capacity	450 lbs.
Maximum Speed	4 mph
Maximum Range	20 miles
Maximum Grade/Incline	8°
Turning Radius	41"
Ground Clearance	4"

Dimensions

Overall Length	48"
Overall Width	26.5" (20" seat) 28.5" (22" seat)
Floor to Seat Height	20"

Seat Dimensions

Available Widths	20" 22"
Depth	20"
Back Height	19"
Armrest to Seat Height	9.5"–13.5"

Power

Motor	Two 24V x 350W x 3800rpm
Controller	Dynamic 50A
Batteries	Two DC12V x 50AH
Battery Charger	11 Amp, Offboard

Brakes

Electromagnetic

Freewheel Mode

Yes

Tires (Flat-free)

Front Casters	9" x 2.5"
Rear Drive Wheels	12.5" x 2.25"

Weights

Without Batteries	100 lbs.
With Batteries	173 lbs.



Feature Focus: Ease of Transport



The Wildcat 450 can be quickly and easily folded to a narrow 13" width, making it perfect for transport. Simply remove the seat cushion and the battery packs.

Accessories

- Elevating Leg Rests
- Swing-away Controller Arm
- Armrest Bag

Color Options

- Black