



## Palm-release Walkers



**ACM7855A**

### Features

- Composite lower side brace stiffens walker frame and adds stability.
- Ergonomic release button and folds easily for transport or storage.
- Folding mechanisms are easy to reach from a seated position.
- Minimal effort is required to open and close the walkers.
- PVC handgrips wear longer and provide greater comfort.
- Available in adult and children's sizes.



### Front Wheel Attachments

- Adding front wheel extensions provides easier maneuverability over surfaces that can be difficult for walker users, while also allowing the user to walk with a more natural stride.
- Available in 5" diameter
- Product Number **ACM7822**.

*Note that wheel attachments are for use on walker front legs ONLY.*

Model	Width Inside Hand Grips	Width Inside Base Legs	Overall Depth	Height Adjustment	Product Weight Capacity
<b>ACM7855A</b>	17"	21.5"	Open: 17.5" Folded: 4.5"	30.375" – 37.375"	300 lbs.
<b>ACM7855C</b>	17"	20.5"	Open: 16.5" Folded: 4.5"	25.375" – 32.375"	300 lbs.

## Dual-release Walkers



**ACM8302**

**ACM8302W**

### Features

- Dual-release folding mechanism allows user to fold walker sides independently
- Height is adjustable using easy snap buttons
- Sturdy, anodized aluminum frame is highly anti-corrosive and lightweight
- Cross brace is located near the hand grips, allowing the user to use a wider, more natural stride
- Foam handgrips for added comfort
- Accommodates all leg attachments
- Rear glide tips allow easier mobility over surfaces that can be difficult for walker users.
- Wheeled version available

Model	Width Inside Hand Grips	Width Inside Base Legs	Overall Depth	Height Adjustment	Product Weight Capacity
<b>ACM8302</b>	18"	22"	17"	32" – 39"	300 lbs.
<b>ACM8302W</b>	18"	22"	17"	32" – 39"	300 lbs.

## Aluminum Bath Bench w/ High Density PE Long-Arc Seat



### Features

- Sturdy, anodized aluminum frame is highly anti-corrosive and lightweight.
- Textured, contoured seat provides comfort and support.
- Plastic seat with molded-in hand holds aid user in rising from seated position.
- Drain holes in seat reduce slipping.
- Angled legs provide additional stability.
- Height is adjustable using easy snap buttons.

Model	Seat Height Range	Seat Depth	Overall Depth	Seat Width	Overall Width	Product Weight Capacity
ACM5101	14" – 22"	12"	17"	20"	20"	300 lbs.

## Aluminum Shower Chair w/ Backrest



### Features

- Sturdy, anodized aluminum frame is highly anti-corrosive and lightweight.
- Textured, molded back and contoured seat provide comfort and support.
- Plastic seat with molded-in hand holds aid user in rising from seated position.
- Drain holes in seat reduce slipping.
- Angled legs for additional stability.
- Easy-to-install back support requires no tools.
- Height is adjustable using easy snap buttons.

Model	Seat Height Range	Seat Depth	Overall Depth	Seat Width	Overall Width	Product Weight Capacity
ACM5104	14" – 22"	12"	17"	20"	20"	300 lbs.

## Aluminum Transfer Bench w/ High Density PE Seat and Backrest



### Features

- Sturdy, anodized aluminum frame is highly anti-corrosive and lightweight.
- Textured, molded back and contoured seat provide comfort and support.
- Seat back with molded-in hand hold for easy carrying.
- Drain holes in seat reduce slipping.
- Angled legs for additional stability.
- Two legs include suction cups for added safety.
- Easy to install back support and armrest tube require no tools.
- Height is adjustable using easy snap buttons.

Model	Seat Height Range	Seat Depth	Overall Depth	Seat Width	Overall Width	Product Weight Capacity
ACM5201N	19.5" – 23.5"	16.75"	20"	27"	30"	300 lbs.

## 3-in-1 Commode



### Features

- Three uses in one convenient unit: Can be used bedside, as toilet safety frame or as a raised toilet seat.
- Seat height easily adjusts from 17" to 25" using snap buttons
- Powder coated steel tube, anti-corrosive construction.
- Textured, blow-molded PE armrest for safety and comfort.
- Non-slip rubber tips for added safety.
- Durable plastic snap-on seat with lid.
- Comes complete with commode bucket, cover and splash shield.

Model	Overall Width	Overall Length	Seat Width	Seat Height Range	Product Weight Capacity
CC1000	25"	16.5"	13.5"	17" – 35"	300 lbs.

## 3-in-1 Foldable Commode



### Features

- Three uses in one convenient unit: Can be used bedside, as toilet safety frame or as a raised toilet seat.
- Seat height easily adjusts from 17" to 25" using snap buttons
- Unit easily disassembles for convenient storage or transport.
- Powder coated steel tube, anti-corrosive construction.
- Textured, blow-molded PE armrest for safety and comfort.
- Non-slip rubber tips for added safety.
- Durable plastic snap-on seat with lid.
- Comes complete with commode bucket, cover and splash shield.

Model	Overall Width	Overall Length	Seat Width	Seat Height Range	Product Weight Capacity
CC2000	25"	16.5"	13.5"	17" – 21"	300 lbs.